

DIAMOND GRIP PLUS™

POWDER-FREE LATEX EXAM GLOVES



• STANDARD EXAMINATION GLOVE



Recommended for use when:

- Durable comfort is required for lasting wear
- Consistent grip is necessary for handling utensils on a day-to-day basis
- Tactile sensitivity is important when working with a number of different instruments



The combination of comfort and strength make **Diamond Grip Plus™** a leading latex glove of choice for many medical professionals. **Diamond Grip Plus™** gloves offer a consistent and secure wet and dry grip. Great elasticity, a comfortable fit, and excellent tactile sensitivity make **Diamond Grip Plus™** a perfect multi-purpose glove.

Diamond Grip Plus™ gloves offer a superior combination of comfort and durability in an all-purpose glove.

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®

REGISTERED
5013485

DIAMOND GRIP PLUS™



POWDER-FREE LATEX EXAM GLOVES



STANDARD EXAMINATION GLOVE

GLOVE PROFILE

	ASTM*	Diamond Grip Plus™
Length (mm/in):	220 / 8.7	248 / 9.8
Cuff Thickness (mils):	n/a	3.9
Palm Thickness (mils):	min. 3.1	5.1
Finger Thickness (mils):	min. 3.1	6.3
Tensile Strength		
Before Aging (MPa):	min. 18	min. 21
After Aging (MPa):	min. 14	min. 18
Elasticity/Elongation:		
Before Aging (%):	min. 650	700
After Aging (%):	min. 500	600
Reduced Protein Rating:	max. 200 µg/dm ²	max. 50 µg/g

GLOVE ATTRIBUTES

Material:	Premium quality natural rubber latex
Type:	Non-sterile
Shape:	Ambidextrous
Color:	Natural
Sizes:	X-small, Small, Medium, Large, X-Large
Cuff:	Beaded
Glove Interior:	Powder-free
Glove Exterior:	Fully textured
Usage:	For single use only
Packaging:	100 gloves per dispenser box, by weight 10 Boxes per case

* ASTM D 3578

Latex Examination Gloves

See for yourself why **Microflex®** gloves are the first choice of exam glove users nationwide. For samples or more information on any of the fine latex and synthetic products from **Microflex®**, contact your **Microflex®** distributor or for the name of a distributor in your area, call:

(800) 876-6866

Or visit our website at:

www.microflex.com

Fax: (800) 876-6632.



PRODUCT CODE	RANGE OF SIZE	DESCRIPTION	CASE QUANTITY
DGP-350-XS	EXTRA-SMALL	Latex, Powder-free, Examination	1000
DGP-350-S	SMALL	Latex, Powder-free, Examination	1000
DGP-350-M	MEDIUM	Latex, Powder-free, Examination	1000
DGP-350-L	LARGE	Latex, Powder-free, Examination	1000
DGP-350-XL	EXTRA-LARGE	Latex, Powder-free, Examination	1000

CAUTION STATEMENT: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
2301 ROBB DRIVE • RENO, NV 89523

©2007 MICROFLEX CORP. ALL RIGHTS RESERVED.
07.05450C AL 060807