

ULTRA ONE®

POWDER-FREE LATEX EXTENDED CUFF EXAMINATION GLOVES

Recommended for use in the following industries:

EMS / Dental / Healthcare / Industrial Safety / Laboratory

For the highest level of protection and durability in latex examination gloves for high-risk environments, rely on **Ultra One**® gloves. Strong, durable and comfortable, **Ultra One**® gloves offer superior defense against many hazards that you come into contact with daily. With an extended cuff for increased protection, the **Ultra One**® glove is a proven performer.

Ultra One® gloves are the best choice in a latex for heavy-duty/high-risk environments.

Recommended for use when:

- Added protection is required for extreme environments
- Extra coverage over the wrist and forearm is needed for enhanced safety and security
- Strength is required for your toughest tasks

MEETS NFPA STANDARDS



Ultra One® Attributes



	ASTM D 3578	Ultra One®
Length (mm/in):	220 / 8.7	300 / 11.8
Cuff Thickness (mils):	n/a	7.1
Palm Thickness (mils):	min. 3.1	9.8
Finger Thickness (mils):	min. 3.1	13.8
Tensile Strength		
Before Aging (MPa):	min. 18	21
After Aging (MPa):	min. 14	18
Elasticity/Elongation:		
Before Aging (%):	min. 650	750
After Aging (%):	min. 500	600
Reduced Protein Rating:	max. 200 µg/dm ²	max. 50 µg/g

Material: Latex	Usage: For single use only
Type: Non-sterile	Shape: Ambidextrous
Color: Natural	Size: S - XL
Cuff: Beaded	Glove Exterior: Textured fingers
Glove Interior: Polymer-coated, powder-free	
Packaging: 500 gloves per case, by weight, 10 boxes per case 50 gloves per dispenser box, by weight	

Order

S ...UL-315-S LUL-315-L
M ...UL-315-M XL ...UL-315-XL

Caution: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. Contact a Microflex® distributor today to purchase quality Microflex® products. For the name of a distributor in your area, call: **(800) 876-6866**. Fax: (800) 876-6632.

