


ColorTouch®

LIGHTLY POWDERED LATEX STANDARD EXAMINATION GLOVES

Recommended for use in the following industries:
Healthcare, Dental



A Colorful Solution in High Quality Hand Protection.

ColorTouch® is a refreshing change from ordinary latex gloves with its peppermint scent, and bright, size-specific colors. Made from premium latex, **ColorTouch®** goes on easily and wears comfortably. Fully textured for a secure grip in wet and dry conditions.

Features and Benefits:

- Premium quality, low modulus latex for greater comfort and reduced hand fatigue
- Lightly powdered with absorbable cornstarch (USP) for easy donning
- Fully textured for a more secure grip
- Meets federal guidelines for maximum powder and maximum water extractable protein content
- Lightly scented with peppermint
- Color-coded sizes

Color Touch's Attributes



	ASTM D 3578	ColorTouch®
Length (mm/in):	220 / 8.7	245 / 9.6
Cuff Thickness (mils):	n/a	3.9
Palm Thickness (mils):	min. 3.1	5.9
Finger Thickness (mils):	min. 3.1	6.7
Tensile Strength		
Before Aging (MPa):	min. 18	21
After Aging (MPa):	min. 14	18
Elasticity/Elongation:		
Before Aging (%):	min. 650	750
After Aging (%):	min. 400	600
Reduced Protein Rating:	max. 200 µg/dm ²	max. 150 µg/g

Order

S . . . CT-133-S Color: Pink
M . . . CT-133-M Color: Green
L CT-133-L Color: Blue

Material: Latex
Type: Non-sterile
Color: Pink, green, blue
Cuff: Beaded
Glove Interior: Lightly powdered with absorbable cornstarch, USP
Packaging: 100 gloves per dispenser box, by weight
10 boxes per case

Usage: For single use only
Shape: Ambidextrous
Size: S - L
Glove Exterior: Fully textured

⚠ Caution: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. Contact a Microflex® distributor today to purchase quality Microflex® products. For the name of a distributor in your area, call: **(800) 876-6866**. Fax: (800) 876-6632.

