

ComfortGrip™

POWDER-FREE LATEX EXAMINATION GLOVES

● STANDARD EXAMINATION GLOVE



Recommended for use when:

- A secure grip is needed in wet and dry environments
- Working with small instruments requires maximum sensitivity
- Comfort, quality, and value are crucial to your needs



If you're looking for the best combination in value and quality, **ComfortGrip™** is the latex glove you need. **ComfortGrip™** gloves are backed by the industry's leading quality system to offer consistent comfort and quality in every box. Fully textured for a secure grip, **ComfortGrip™** gloves provide the fit and protection you need at a competitive price.

ComfortGrip™ gloves offer the best overall value and quality in a latex glove.

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®

REGISTERED
ISO 13485

ComfortGrip™

POWDER-FREE LATEX EXAMINATION GLOVES

STANDARD EXAMINATION GLOVE



GLOVE PROFILE

	ASTM*	ComfortGrip™
Length (mm/in):	220 / 8.7	245 / 9.6
Cuff Thickness (mils):	n/a	3.9
Palm Thickness (mils):	min. 3.1	5.1
Finger Thickness (mils):	min. 3.1	5.9
Tensile Strength		
Before Aging (MPa):	min. 18	min. 18
After Aging (MPa):	min. 14	min. 14
Elasticity/Elongation:		
Before Aging (%):	min. 650	650
After Aging (%):	min. 500	500
Reduced Protein Rating:	max. 200 µg/dm ²	max. 50 µg/g

GLOVE ATTRIBUTES

Material:	Natural rubber latex
Type:	Non-sterile
Shape:	Ambidextrous
Color:	Natural
Sizes:	X-small, Small, Medium, Large, X-Large
Cuff:	Beaded
Glove Interior:	Powder-free, chlorinated, polymer coated
Glove Exterior:	Fully textured
Usage:	For single use only
Packaging:	100 gloves per dispenser box, by weight 10 Boxes per case

* ASTM D 3578

Latex Examination Gloves

See for yourself why **Microflex®** gloves are the first choice of exam glove users nationwide. For samples or more information on any of the fine latex and synthetic products from **Microflex®**, contact your **Microflex®** distributor or for the name of a distributor in your area, call:

(800) 876-6866

Or visit our website at:

www.microflex.com

Fax: (800) 876-6632.



PRODUCT CODE	RANGE OF SIZE	DESCRIPTION	CASE QUANTITY
CFG-900-XS	EXTRA-SMALL	Latex, Powder-free, Examination	1000
CFG-900-S	SMALL	Latex, Powder-free, Examination	1000
CFG-900-M	MEDIUM	Latex, Powder-free, Examination	1000
CFG-900-L	LARGE	Latex, Powder-free, Examination	1000
CFG-900-XL	EXTRA-LARGE	Latex, Powder-free, Examination	1000

CAUTION STATEMENT: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.