

Microflex®
SUPRENO® SE

Powder-Free Nitrile
 Exam Gloves

Item #

Extra Small.....	SU-690-XS
Small.....	SU-690-S
Medium.....	SU-690-M
Large.....	SU-690-L
Extra Large.....	SU-690-XL



Recommended for use in the following industries:

Automotive / Dental / EMS / Healthcare / Industrial Safety / Laboratory

If your need is strength and durability, **Supreno® SE** gloves provide exceptional toughness in a comfortable nitrile glove. **Supreno® SE** gloves provide reliable protection against the hazards you work with day in and day out.

Supreno® SE gloves are the best combination of strength and durability in an all-purpose synthetic glove.

Supreno® SE gloves meet or exceed the emergency medical glove requirements of NFPA 1999, Standard on Protective Clothing for Emergency Medical Operations, 2008 edition.¹

¹Certification on file.

Recommended for use when:

- Strength and protection are of great importance
- Day-to-day hazards require the best combination of durability and security
- Working in environments that demand consistent grip and reliability

Microflex® Supreno® SE Profile

For current specification information, please visit www.microflex.com.

	ASTM D 6319	Supreno® SE		
Length (mm/in):	230 / 9.1	245 / 9.6	Material: Nitrile	Type: Non-sterile
Cuff Thickness (mils):	n/a	3.9	Size: XS - XL	Color: Blue-Violet
Palm Thickness (mils):	min. 2.0	4.7	Shape: Ambidextrous	Usage: For single use only
Finger Thickness (mils):	min. 2.0	7.1	Cuff: Standard examination length (9.5"), Beaded	
Tensile Strength:			Glove Exterior: Textured Fingers	
Before Aging (MPa):	min. 14	30	Glove Interior: Powder-Free, Chlorinated	
After Aging (MPa):	min. 14	30	Packaging: 1,000 gloves/case, by weight; 10 boxes/case	
Elasticity/Elongation:			100 gloves per dispenser box, by weight	
Before Aging (%):	min. 500	500		
After Aging (%):	min. 400	400	Caution: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.	

The reason is the quality... the difference is Microflex®.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit www.microflex.com. For the name of a distributor in your area, call: **(800) 876-6866**.

