

POWERGRIP®


LIGHTLY POWDERED LATEX EXAMINATION GLOVES

• STANDARD EXAMINATION GLOVE



Features & Benefits:

- Low modulus latex glove allows for reduced hand fatigue over long periods of wear
- Fully textured for a secure grip in wet or dry conditions
- Longer than typical latex examination gloves
- Lightly powdered with absorbable cornstarch (USP) for easy donning



One Powerful Performer.

Lightly powdered to go on easily, this sturdy examination glove is fully textured for professionals who need a firm grip in wet or dry conditions. Greater thickness and a longer cuff offer coverage over the wrist and forearm, providing increased protection where it is most needed. **PowerGrip®** meets or exceeds proposed FDA and current ASTM standards for exam glove quality and performance.

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®

REGISTERED
ISO13485

POWERGRIP®

LIGHTLY POWDERED LATEX EXAMINATION GLOVES

STANDARD EXAMINATION GLOVE

GLOVE PROFILE

	ASTM*	PowerGrip®**
Length (mm/in):	220 (8.7) / 230 (9.1)	240/9.6
Cuff Thickness (mils):	n/a	3.9
Palm Thickness (mils):	min. 3.1	5.9
Finger Thickness (mils):	min. 3.1	6.7
Tensile Strength		
Before Aging (MPa):	min. 18	21
After Aging (MPa):	min. 14	18
Elasticity/Elongation:		
Before Aging (%):	min. 650	750
After Aging (%):	min. 500	600
Reduced Protein Rating:	max. 200 µg/dm ²	max. 150 µg/g

GLOVE ATTRIBUTES

Material:	Premium quality natural rubber latex
Type:	Non-sterile
Shape:	Ambidextrous
Color:	White
Sizes:	X-Small, Small, Medium, Large, X-Large
Cuff:	Beaded
Glove Interior:	Lightly powdered with absorbable cornstarch, USP
Glove Exterior:	Textured
Usage:	For single use only
Packaging:	100 gloves per dispenser box, by weight 10 Boxes per case

* ASTM D 3578



Latex Examination Gloves

See for yourself why **Microflex®** gloves are the first choice of exam glove users nationwide. For samples or more information on any of the fine latex and synthetic products from **Microflex®**, contact your **Microflex®** distributor or for the name of a distributor in your area, call:

(800) 876-6866

Or visit our website at:

www.microflex.com

Fax: (800) 876-6632.

PRODUCT CODE	RANGE OF SIZE	DESCRIPTION	CASE QUANTITY
PG-199-XS	EXTRA-SMALL	Latex, Lightly Powdered, Examination	1000
PG-199-S	SMALL	Latex, Lightly Powdered, Examination	1000
PG-199-M	MEDIUM	Latex, Lightly Powdered, Examination	1000
PG-199-L	LARGE	Latex, Lightly Powdered, Examination	1000
PG-199-XL	EXTRA-LARGE	Latex, Lightly Powdered, Examination	1000

CAUTION STATEMENT: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.



MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
2301 ROBB DRIVE • RENO, NV 89523