

Microflex®
ColorTouch®

Lightly Powdered
Latex Exam Gloves

Item #

Small CT-133-S
 Medium..... CT-133-M
 Large..... CT-133-L



Recommended for use in the following industries:

Dental / Healthcare / Industrial Safety / Laboratory / Tattoo

A Colorful Solution in High Quality Hand Protection.

ColorTouch® gloves are a refreshing change from ordinary latex gloves, with a peppermint scent and bright, size-specific colors. Made from premium latex, **ColorTouch®** gloves go on easily and wear comfortably. Fully textured for a secure grip in wet and dry conditions.

- Lightly powdered with absorbable cornstarch (USP) for easy donning
- Fully textured for a more secure grip
- Lightly scented with peppermint
- Color-coded sizes

Microflex® ColorTouch® Profile

For current specification information, please visit www.microflex.com.

	ASTM D 3578	ColorTouch®		
Length (mm/in):	230 / 9.1	245 / 9.6	Material: Latex	Type: Non-sterile
Cuff Thickness (mils):	n/a	3.9	Size: XS - L	Color: Pink, Green, Blue
Palm Thickness (mils):	min. 3.1	5.9	Shape: Ambidextrous	Usage: For single use only
Finger Thickness (mils):	min. 3.1	6.7	Cuff: Standard examination length (9.5"), Beaded	
Tensile Strength:			Glove Exterior: Fully Textured	
Before Aging (MPa):	min. 18	22	Glove Interior: Lightly Powdered, absorbable cornstarch (USP)	
After Aging (MPa):	min. 14	21	Packaging: 1,000 gloves/case, by weight; 10 boxes/case	
Elasticity/Elongation:			100 gloves per dispenser box, by weight	
Before Aging (%):	min. 650	750		
After Aging (%):	min. 500	600	Caution: This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.	

The reason is the quality... the difference is Microflex®.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit www.microflex.com. For the name of a distributor in your area, call: **(800) 876-6866**.

MICROFLEX®
 THE MOST TRUSTED NAME IN GLOVES®
 A division of BarrierSafe Solutions International, Inc.