

Microflex®

Derma Free®

Powder-Free Vinyl Exam Gloves

Item

Small DF-850-S
 Medium..... DF-850-M
 Large..... DF-850-L
 Extra Large..... DF-850-XL



Recommended for use in the following industries:

Automotive / Dental / EMS / Healthcare / Industrial Safety / Laboratory

Comfort, performance and value in a powder-free medical examination glove.

Backed by the industry's leading quality system, **Derma Free®** powder-free vinyl gloves provide solid hand protection for a variety of applications. Thicker than typical vinyl gloves, the **Derma Free®** glove is 100% synthetic and features a smooth surface for more natural tactile sensitivity.

Backed by the industry's leading quality system, **Derma Free®** powder-free vinyl is 100% synthetic and features a smooth surface for natural tactile sensitivity.

Microflex® Derma Free® Profile

For current specification information, please visit www.microflex.com.

	ASTM D 5250	Derma Free®		
Length (mm/in):	230 / 9.1	230 / 9.1	Material: Vinyl	Type: Non-sterile
Cuff Thickness (mils):	n/a	min. 2.0	Size: XS - XL	Color: Clear
Palm Thickness (mils):	min. 3.1	min. 3.1	Shape: Ambidextrous	Usage: For single use only
Finger Thickness (mils):	min. 2.0	min. 2.0	Cuff: Standard examination length (9"), Beaded	
Tensile Strength:			Glove Exterior: Smooth	
Before Aging (MPa):	min. 11	13	Glove Interior: Powder-Free	
After Aging (MPa):	min. 11	13	Packaging: 1,000 gloves/case, by weight; 10 boxes/case	
Elasticity/Elongation:			100 gloves per dispenser box, by weight	
Before Aging (%):	min. 300	300		
After Aging (%):	min. 300	300	Caution: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.	

The reason is the quality... the difference is Microflex®.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit www.microflex.com. For the name of a distributor in your area, call: **(800) 876-6866**.

MICROFLEX®
 THE MOST TRUSTED NAME IN GLOVES®
 A division of BarrierSafe Solutions International, Inc.