

Microflex®

# EVOLUTION ONE®

## Powder-Free Latex Exam Gloves

### Item #

Extra Small.....	EV-2050-XS
Small.....	EV-2050-S
Medium.....	EV-2050-M
Large.....	EV-2050-L
Extra Large.....	EV-2050-XL



### Recommended for use in the following industries:

Dental / Healthcare / Industrial Safety / Laboratory

Premium comfort and consistency. **Evolution One®** gloves offer a secure grip and don easily in wet and dry environments.

**Evolution One®** gloves provide maximum comfort combined with enhanced tactile sensitivity, durability and strength.

**Evolution One®** gloves are made with a premium latex for maximum comfort and consistency in mid-range environments.

### Recommended for use when:

- Maximum comfort is required for lasting wear
- Rapid donning and anti-sticking properties are needed
- Working with a variety of instruments and equipment

### Microflex® Evolution One® Profile

For current specification information, please visit [www.microflex.com](http://www.microflex.com).

	ASTM D 3578	Evolution One®		
<b>Length (mm/in):</b>	230 / 9.1	245 / 9.6	<b>Material:</b> Latex	<b>Type:</b> Non-sterile
<b>Cuff Thickness (mils):</b>	n/a	3.9	<b>Size:</b> XS - XL	<b>Color:</b> Natural
<b>Palm Thickness (mils):</b>	min. 3.1	5.5	<b>Shape:</b> Ambidextrous	<b>Usage:</b> For single use only
<b>Finger Thickness (mils):</b>	min. 3.1	5.9	<b>Cuff:</b> Standard examination length (9.5"), Beaded	
<b>Tensile Strength:</b>			<b>Glove Exterior:</b> Fully Textured, Polymer-Coated	
Before Aging (MPa):	min. 18	24	<b>Glove Interior:</b> Powder-Free, Chlorinated	
After Aging (MPa):	min. 14	18	<b>Packaging:</b> 1,000 gloves/case, by weight; 10 boxes/case	
<b>Elasticity/Elongation:</b>			100 gloves per dispenser box, by weight	
Before Aging (%):	min. 650	700		
After Aging (%):	min. 500	600		

**Caution:** This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

### The reason is the quality... the difference is Microflex®.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit [www.microflex.com](http://www.microflex.com). For the name of a distributor in your area, call: **(800) 876-6866**.

**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®  
A division of BarrierSafe Solutions International, Inc.