

Microflex®
FreeForm® EC

Powder-Free Nitrile
 Extended Cuff Exam Gloves

Item #

- Small FFE-775-S
- Medium..... FFE-775-M
- Large..... FFE-775-L
- Extra Large..... FFE-775-XL
- XX-Large..... FFE-775-XXL



Recommended for use in the following industries:

EMS / Healthcare / Industrial Safety / Laboratory

With an extended cuff for added protection over the wrist and arm, **FreeForm® EC** gloves offer comfortable strength in a relaxed fit nitrile examination product.

Unlike conventional nitrile gloves, which tend to feel stiff and unforgiving, **FreeForm® EC** gloves work with your hands as they move. **FreeForm® EC** gloves are a model example of strength and comfort in a nitrile product.

FreeForm® EC gloves meet or exceed the emergency medical glove requirements of NFPA 1999, Standard on Protective Clothing for Emergency Medical Operations, 2008 edition.*

Recommended for use:

- Strength and protection are of great importance
- Day-to-day hazards require the best combination of durability and comfort
- Working in environments that demand consistent grip and reliability
- Non-stick properties that are resistant to tapes and adhesives

*Certification on file.

Microflex® FreeForm® EC Profile			For current specification information, please visit www.microflex.com .	
	ASTM D 6319	FreeForm® EC		
Length (mm/in):	220 / 8.7	290 / 11.4	Material: Nitrile	Type: Non-sterile
Cuff Thickness (mils):	n/a	3.5	Size: S - XXL	Color: Blue
Palm Thickness (mils):	min. 2.0	4.7	Shape: Ambidextrous	Usage: For single use only
Finger Thickness (mils):	min. 2.0	6.3	Cuff: Extended Cuff length (11.4"), Beaded	
Tensile Strength:			Glove Exterior: Textured Fingers	
Before Aging (MPa):	min. 14	31	Glove Interior: Powder-Free, Chlorinated	
After Aging (MPa):	min. 14	27	Packaging: 500 gloves/case, by weight; 10 boxes/case	
Elasticity/Elongation:			50 gloves per dispenser box, by weight	
Before Aging (%):	min. 500	500	Caution: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.	
After Aging (%):	min. 400	400		

The reason is the quality... the difference is Microflex®.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit www.microflex.com. For the name of a distributor in your area, call: **(800) 876-6866**.

