

Microflex®  
**MICRO ONE®**

Lightly Powdered  
**Latex Exam Gloves**

**Item #**

Extra Small.....	MO-150-XS
Small .....	MO-150-S
Medium.....	MO-150-M
Large.....	MO-150-L
Extra Large.....	MO-150-XL



**Recommended for use in the following industries:**

Automotive / Dental / Healthcare / Industrial Safety / Laboratory / Tattoo

**Lightly powdered. A Proven Performer.**

**Micro One®** gloves are a highly versatile examination glove that performs well in demanding conditions. Made from soft latex, **Micro One®** gloves have a natural fit for reduced hand fatigue. Lightly powdered **Micro One®** gloves don easily and have a smooth surface for uncompromising sensitivity. **Micro One®** gloves meet or exceed current ASTM standards for exam glove quality and performance.

**Recommended for use when:**

- Lightly powdered with absorbable cornstarch (USP) for easy donning
- Smooth surface for excellent tactile sensitivity
- Meets federal guidelines for maximum powder and maximum water extractable protein content

**Microflex® Micro One® Profile**

For current specification information, please visit [www.microflex.com](http://www.microflex.com).

	ASTM D 3578	Micro One®		
<b>Length (mm/in):</b>	230 / 9.1	245 / 9.6	<b>Material:</b> Latex	<b>Type:</b> Non-sterile
<b>Cuff Thickness (mils):</b>	n/a	3.9	<b>Size:</b> XS - XL	<b>Color:</b> White
<b>Palm Thickness (mils):</b>	min. 3.1	5.5	<b>Shape:</b> Ambidextrous	<b>Usage:</b> For single use only
<b>Finger Thickness (mils):</b>	min. 3.1	6.3	<b>Cuff:</b> Standard examination length (9.5"), Beaded	
<b>Tensile Strength:</b>			<b>Glove Exterior:</b> Smooth	
Before Aging (MPa):	min. 18	27	<b>Glove Interior:</b> Lightly Powdered, absorbable cornstarch (USP)	
After Aging (MPa):	min. 14	24	<b>Packaging:</b> 1,000 gloves/case, by weight; 10 boxes/case	
<b>Elasticity/Elongation:</b>			100 gloves per dispenser box, by weight	
Before Aging (%):	min. 650	750	<b>Caution:</b> This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.	
After Aging (%):	min. 500	600		

**The reason is the quality... the difference is Microflex®.**

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit [www.microflex.com](http://www.microflex.com). For the name of a distributor in your area, call: **(800) 876-6866**.

