


Microflex®  
**SUPRENO® EC**

Powder-Free Nitrile  
 Extended Cuff Exam Gloves

**Item #**

Small .....	SEC-375-S
Medium .....	SEC-375-M
Large .....	SEC-375-L
Extra Large .....	SEC-375-XL
XX-Large .....	SEC-375-XXL
XXX-Large .....	SEC-375-XXXL



**Recommended for use in the following industries:**

Dental / EMS / Funeral / Healthcare / Industrial Safety / Laboratory

**Exceptional Protection in a Tough Nitrile Glove.**

**Supreno® EC** gloves are a durable nitrile glove for high-risk applications. **Supreno® EC** gloves offer all the benefits of a nitrile glove, with the most important being unique and reliable protection against the hazards you work with. **Supreno® EC** gloves are formulated for enhanced strength and durability while maintaining a comfortable fit.

**Supreno® EC** gloves meet or exceed the emergency medical glove requirements of NFPA 1999, Standard on Protective Clothing for Emergency Medical Operations, 2008 edition.<sup>1</sup>

*Tested for use with chemotherapy drugs. Gloves used for protection against chemotherapy drugs exposure must be selected specifically for the type of chemicals used. Tested for use against Vincristine, Dacarbazine, Amethopterin Hydrate, Carmustine, Cyclophosphamide, Doxorubicin Hydrochloride, 5-Fluorouracil, Cisplatin, Etoposide, Paclitaxel, Thio-Tepa.*

<sup>1</sup>Certification on file.

**Recommended for use:**

- Strength and protection are of great importance
- Extra coverage over the wrist and forearm is needed for enhanced safety and security
- Day-to-day hazards require the best combination of durability and security
- Working in environments that demand consistent grip and reliability
- Non-stick properties that are resistant to tapes and adhesives

**Microflex® Supreno® EC Profile**

For current specification information, please visit [www.microflex.com](http://www.microflex.com).

	ASTM D 6319	Supreno® EC		
<b>Length (mm/in):</b>	220 / 8.7	295 / 11.6	<b>Material:</b> Nitrile	<b>Type:</b> Non-sterile
<b>Cuff Thickness (mils):</b>	n/a	4.3	<b>Size:</b> S - XXXL	<b>Color:</b> Blue-Violet
<b>Palm Thickness (mils):</b>	min. 2.0	5.5	<b>Shape:</b> Ambidextrous	<b>Usage:</b> For single use only
<b>Finger Thickness (mils):</b>	min. 2.0	7.9	<b>Cuff:</b> Extended Cuff length (12"), Beaded	
<b>Tensile Strength:</b>			<b>Glove Exterior:</b> Textured Fingers, Polymer-Coated	
Before Aging (MPa):	min. 14	30	<b>Glove Interior:</b> Powder-Free, Chlorinated	
After Aging (MPa):	min. 14	27	<b>Packaging:</b> 500 gloves/case, by weight; 10 boxes/case	
<b>Elasticity/Elongation:</b>			50 gloves per dispenser box, by weight	
Before Aging (%):	min. 500	500	XXX-Large: 40 gloves per dispenser box, by weight	
After Aging (%):	min. 400	400	<b>Caution:</b> Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.	

**The reason is the quality... the difference is Microflex®.**

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit [www.microflex.com](http://www.microflex.com). For the name of a distributor in your area, call: **(800) 876-6866**.

