

Microflex®

# SYNETRON®

## Powder-Free Latex Extended Cuff Exam Gloves

### Item #

- Small ..... SY-911-S
- Medium ..... SY-911-M
- Large ..... SY-911-L
- Extra Large ..... SY-911-XL
- XX-Large ..... SY-911-XXL
- XXX-Large ..... SY-911-XXXL



### Recommended for use in the following industries:

EMS / Industrial Safety / Laboratory / Law & Corrections

Textured for a secure grip. **Synetron®** gloves are polymer-coated for easy donning. This glove meets the single-use emergency medical examination glove requirements of NFPA 1999, Standard on Protective Clothing for Emergency Medical Operations, 2008 edition.\*

- Interior/exterior polymer system provides superior anti-sticking properties resulting in highly effective donning in wet, damp or dry conditions
- Extended cuff for increased protection over the wrist and forearm
- Fine fingertip sensitivity for finding a pulse and starting an I.V.

\*Certification on file.

Microflex® Synetron® Profile			For current specification information, please visit <a href="http://www.microflex.com">www.microflex.com</a> .	
	ASTM D 3578	Synetron®		
<b>Length (mm/in):</b>	230 / 9.1	290 / 11.4	<b>Material:</b> Latex	<b>Type:</b> Non-sterile
<b>Cuff Thickness (mils):</b>	n/a	6.7	<b>Size:</b> S - XXXL	<b>Color:</b> Natural
<b>Palm Thickness (mils):</b>	min. 3.1	9.1	<b>Shape:</b> Ambidextrous	<b>Usage:</b> For single use only
<b>Finger Thickness (mils):</b>	min. 3.1	11.8	<b>Cuff:</b> Extended Cuff length (12"), Beaded	
<b>Tensile Strength:</b>			<b>Glove Exterior:</b> Fully Textured, Polymer-Coated	
Before Aging (MPa):	min. 18	26	<b>Glove Interior:</b> Powder-Free, Chlorinated	
After Aging (MPa):	min. 14	22	<b>Packaging:</b> 500 gloves/case, by weight; 10 boxes/case	
<b>Elasticity/Elongation:</b>			50 gloves per dispenser box, by weight	
Before Aging (%):	min. 650	700	XXX-Large: 40 gloves per dispenser box, by weight	
After Aging (%):	min. 500	600	<b>Caution:</b> This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.	

### The reason is the quality... the difference is Microflex®.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit [www.microflex.com](http://www.microflex.com). For the name of a distributor in your area, call: **(800) 876-6866**.

**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®  
A division of BarrierSafe Solutions International, Inc.