

Microflex®  
**ULTRA ONE®**  
 Powder-Free Latex  
 Extended Cuff Exam Gloves

**Item #**  
 Small ..... UL-315-S  
 Medium..... UL-315-M  
 Large..... UL-315-L  
 Extra Large..... UL-315-XL



**Recommended for use in the following industries:**  
 Automotive / EMS / Industrial Safety / Laboratory

For the highest level of protection and durability in latex examination gloves for high-risk environments, rely on **Ultra One®** gloves. Strong, durable and comfortable, **Ultra One®** gloves offer superior defense against many hazards that you come into contact with daily. With an extended cuff for increased protection, the **Ultra One®** glove is a proven performer.

**Ultra One®** gloves are the best choice in a latex for heavy-duty/high-risk environments.

- Recommended for use when:**
- Added protection is required for extreme environments
  - Extra coverage over the wrist and forearm is needed for enhanced safety and security
  - Strength is required for your toughest tasks

NFPA certification on file.

Microflex® Ultra One® Profile			For current specification information, please visit <a href="http://www.microflex.com">www.microflex.com</a> .	
	ASTM D 3578	Ultra One®	<b>Material:</b> Latex	<b>Type:</b> Non-sterile
<b>Length (mm/in):</b>	230 / 9.1	300 / 11.8	<b>Size:</b> S - XL	<b>Color:</b> Natural
<b>Cuff Thickness (mils):</b>	n/a	6.7	<b>Shape:</b> Ambidextrous	<b>Usage:</b> For single use only
<b>Palm Thickness (mils):</b>	min. 3.1	9.8	<b>Cuff:</b> Extended Cuff length (12"), Beaded	
<b>Finger Thickness (mils):</b>	min. 3.1	11.8	<b>Glove Exterior:</b> Textured Fingers	
<b>Tensile Strength:</b>			<b>Glove Interior:</b> Powder-Free, Chlorinated	
Before Aging (MPa):	min. 18	33	<b>Packaging:</b> 500 gloves/case, by weight; 10 boxes/case	
After Aging (MPa):	min. 14	25	50 gloves per dispenser box, by weight	
<b>Elasticity/Elongation:</b>			<b>Caution:</b> This product contains natural rubber latex (latex) which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.	
Before Aging (%):	min. 650	750		
After Aging (%):	min. 500	600		

**The reason is the quality... the difference is Microflex®.**  
 See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit [www.microflex.com](http://www.microflex.com). For the name of a distributor in your area, call: **(800) 876-6866**.

