

Microflex®

UltraSense® Powder-Free Blue Nitrile Exam Gloves

Item

Extra Small.....	US-220-XS
Small	US-220-S
Medium.....	US-220-M
Large.....	US-220-L
Extra Large.....	US-220-XL



Recommended for use in the following industries:

Dental / EMS / Healthcare / Industrial Safety / Laboratory

A nitrile glove with the fit and feel of latex.

UltraSense® gloves are made from a next-generation nitrile formulation to provide comfort, while maintaining the reliable barrier protection of a synthetic glove. **UltraSense®** gloves provides a precise grip in both wet and dry conditions, while offering comfort for extended periods of use.

Recommended for use:

- Handling small instruments, where fit and feel is key
- Light tasks require excellent grip on a day-to-day basis
- Performing extended or repetitive duties that may cause hand fatigue
- Non-stick properties that are resistant to tapes and adhesives

Microflex® UltraSense® Profile

For current specification information, please visit www.microflex.com.

	ASTM D 6319	UltraSense®		
Length (mm/in):	230 / 9.1	243 / 9.6	Material: Nitrile	Type: Non-sterile
Cuff Thickness (mils):	n/a	2.8	Size: XS - XL	Color: Blue
Palm Thickness (mils):	min. 2.0	3.1	Shape: Ambidextrous	Usage: For single use only
Finger Thickness (mils):	min. 2.0	4.3	Cuff: Standard examination length (9.5"), Beaded	
Tensile Strength:			Glove Exterior: Textured Fingers, Polymer-Coated	
Before Aging (MPa):	min. 14	30	Glove Interior: Powder-Free, Chlorinated	
After Aging (MPa):	min. 14	25	Packaging: 1,000 gloves/case, by weight; 10 boxes/case	
Elasticity/Elongation:			100 gloves per dispenser box, by weight	
Before Aging (%):	min. 500	500	Caution: Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.	
After Aging (%):	min. 400	400		

The reason is the quality... the difference is Microflex®.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit www.microflex.com. For the name of a distributor in your area, call: **(800) 876-6866**.

MICROFLEX®
THE MOST TRUSTED NAME IN GLOVES®
A division of BarrierSafe Solutions International, Inc.