

Microflex®

# UltraSense® EC

## Powder-Free Nitrile Extended Cuff Exam Gloves

### Item #

- Extra Small..... USE-880-XS
- Small ..... USE-880-S
- Medium..... USE-880-M
- Large..... USE-880-L
- Extra Large..... USE-880-XL
- XX-Large..... USE-880-XXL
- XXX-Large..... USE-880-XXXL



### Recommended for use in the following industries:

Dental / EMS / Healthcare / Industrial Safety / Laboratory

### The protection you deserve at a price you can afford!

**UltraSense® EC** gloves are a next generation nitrile formulation that provides a reliable barrier while offering the comfort you are used to with a latex glove. Engineered with enhanced tactile sensitivity that provides an excellent grip in both wet and dry environments. **UltraSense® EC** gloves are an extended cuff offering that is packed 100 gloves per box as opposed to the standard 50, maximizing your savings.

### Recommended for use:

- Strength and protection are of great importance
- Working in environments that demand consistent grip and reliability
- Precise contact requires enhanced tactile sensitivity
- Non-stick properties that are resistant to tapes and adhesives

### Microflex® UltraSense® EC Profile

For current specification information, please visit [www.microflex.com](http://www.microflex.com).

	ASTM D 6319	UltraSense® EC		
<b>Length (mm/in):</b>	230 / 9.1	290 / 11.4	<b>Material:</b> Nitrile	<b>Type:</b> Non-sterile
<b>Cuff Thickness (mils):</b>	n/a	2.8	<b>Size:</b> XS - XXXL	<b>Color:</b> Blue
<b>Palm Thickness (mils):</b>	min. 2.0	3.1	<b>Shape:</b> Ambidextrous	<b>Usage:</b> For single use only
<b>Finger Thickness (mils):</b>	min. 2.0	4.3	<b>Cuff:</b> Extended Cuff length (12"), Beaded	
<b>Tensile Strength:</b>			<b>Glove Exterior:</b> Textured Fingers, Polymer-Coated	
Before Aging (MPa):	min. 14	33	<b>Glove Interior:</b> Powder-Free, Chlorinated	
After Aging (MPa):	min. 14	32	<b>Packaging:</b> 1,000 gloves/case, by weight; 10 boxes/case	
<b>Elasticity/Elongation:</b>			100 gloves per dispenser box, by weight	
Before Aging (%):	min. 500	500		
After Aging (%):	min. 400	400		

**Caution:** Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for use.

### The reason is the quality... the difference is Microflex®.

See for yourself why Microflex® gloves are the first choice of exam glove users nationwide. For **FREE SAMPLES**, please visit [www.microflex.com](http://www.microflex.com). For the name of a distributor in your area, call: **(800) 876-6866**.

**MICROFLEX®**  
THE MOST TRUSTED NAME IN GLOVES®  
A subsidiary of BarrierSafe Solutions International, Inc.